

Pin Liner Instructions



First completely turn the insert inside out.



Place the gel against the bottom of the residual limb making sure to place the pin in the center of the limb with the white side or logo on the front side.



Begin to roll up the limb making sure to avoid wrinkles as you roll up the liner.



Make sure the white portion of the liner is on the front of the leg as shown.

Pin placement should be in the center of the limb.

Pin Liner Instructions



When using an insert in conjunction to the roll on liner, place it on over the gel liner as shown with the pin through the distal end of the insert.

Distal view of the insert. Pin placement center hole.



Once you have the insert in its proper placement, place the insert into the laminated socket and push with firm pressure until you hear at least one click noise. This sound means the lock has secured the liner into the socket. Now you're ready to stand and walk safely.



To remove the limb from the socket you must place your thumb and hold pressure on the release pin as shown. You must keep firm pressure until the limb has fully been removed from the socket.